

IMPLEMENTATION OF ECG AND PULSE-OXIMETER

- Under the guidance of
MR.K.DAYANAND,M.TECH
ASSISTANT PROFESSOR

- Presented by
T.LAHARI-20FE1A04H1
S.AMEENA-20FE1A04E1
P.AKASH-20FE1A04C6
T.VAMSI-20FE1A04H0



VIGNAN'S LARA
INSTITUTE OF TECHNOLOGY & SCIENCE

Approved by AICTE New Delhi & Affiliated to JNTUK Kakinada

An ISO 9001 : 2015 Certified Institution

CONTENTS

- ABSTRACT
- INTRODUCTION
- COMPONENTS REQUIRED
- WORKING PRINCIPLE
- SOURCE CODE
- ADVANTAGES AND APPLICATIONS
- ADVANCEMENTS
- EXPECTED RESULT
- CONCLUSION

ABSTRACT

- The aim of the ECG is to detect the abnormal heart rhythms
- An ECG is to help diagnose and monitor conditions affecting the heart. It can be used to investigate symptoms of a possible heart problem, such as chest pain, palpitations, dizziness and shortness of breath.
- Pulse oximetry is a test used to measure the oxygen level (oxygen saturation) of the blood. It is an easy, painless measure of how well oxygen is being sent to parts of your body furthest from your heart, such as the arms and legs.


INTRODUCTION

➤ PULSE OXIMETER:



- It is a non-invasive device used to measure blood oxygen levels and can also display heartbeat
- A pulse oximeter is a MEDICAL DEVICE that indirectly monitors the oxygen saturation of a patient's blood
- Pulse oximeter measures:
 1. The oxygen saturation of haemoglobin in arterial blood
 2. The pulse rate in beats per minute
- It is a safe method and a simple method of assessing oxygenation
- Convenient and measurement can be continuous
- The maintenance of optimal O₂ delivery is the core concern during anaesthesia

ELECTROCARDIOGRAPH:



shutterstock.com • 246876961


- The ecg machine is designed to recognise and record any electrical activity within the heart
- It provide information about the function of the intracardiac conducting tissue of the heart and reflects the presence of cardiac disease through its electrical properties
- With each heartbeat, an electrical impulse starts from the superior part of the heart to the bottom. the impulse prompts the heart to contract and pumps blood
- Understanding ecg helps to understand how the heart works
- Electrical activity has two phases:
 1. phase of depolarisation
 2. phase of repolarisation

COMPONENTS REQUIRED

- Aurdino uno
- AD8232 ECG SENSOR
- MX30100 SENSOR
- NODEMCU ESP8266
- JUMPER WIRES


ARDUINO UNO:

- The arduino uno is an open source microcontroller board based on the microchip ATmega328P microcontroller and developed by aurdino.cc and initially realeased in 2010.
- The board is equipped with sets of digital and analog inputur output pins that may be interfaced to various expansion board and circuits.



AD8232 ECG SENSOR:

- The AD8232 ECG sensor is a commercial board used to calculate the electrical movement of the human heart.
- Electrocardiograms can be very noisy, so to reduce the noise the AD8232 chip can be used.
- The working principle of the ECG sensor is like an operational amplifier to help in getting a clear signal from the intervals simply.




MAX30100 SENSOR:

- MAX30100 is an integrated pulse oximeter and heart rate monitor sensor solution.
- It combines two LED's photodetector, optimized optics, and low-noise analog signal processing to detect pulse oximeter and heart-rate signals.
- The sensor has two light-emitting diodes and one photodiode. The LED's are used to emit the light and the photodiode is used to detect and measure the intensity of the received light.




NODEMCU:

- NodeMCU is an open source platform based on ESP8266 which can connect objects and let data transfer using the Wi-Fi protocol.
- In addition, by providing some of the most important features of microcontrollers such as GPIO, PWM, ADC .
- The advantages of NodeMCU are low cost, integrated support for WiFi networks, a smaller board size, and lower energy




WORKING PRINCIPLE OF ECG:



- It works on the principle that a contracting muscle generates a small electric current that can be detected and measured through electrodes suitably placed on the body.
- The electrode picks up the current and transmit them to an amplifier inside the electrocardiograph. Then electrocardiograph amplifies the current and records them on a paper as a wavy line.
- In an electrocardiograph, a sensitive lever traces the changes in current on a moving sheet of paper.
- A modern electrocardiograph may also be connected to an oscilloscope, an instrument that display the current on a screen.

WORKING PRINCIPLE OF PULSE OXIMETER:



- The device has two LEDs, one emitting red light, another emitting infrared light. For pulse rate, only the infrared light is needed. Both the red light and infrared light is used to measure oxygen levels in the blood.
- When the heart pumps blood, there is an increase in oxygenated blood as a result of having more blood. As the heart relaxes, the volume of oxygenated blood also decreases. By knowing the time between the increase and decrease of oxygenated blood, the pulse rate is determined.
- It turns out, oxygenated blood absorbs more infrared light and passes more red light while deoxygenated blood absorbs red light and passes more infrared light.

SOURCE CODE OF ECG

- `void setup() {`
- `// initialize the serial communication:`
- `Serial.begin(9600);`
- `pinMode(10, INPUT); // Setup for leads off detection LO +`
- `pinMode(11, INPUT); // Setup for leads off detection LO -`
- `}`
- `void loop() {`
- `if((digitalRead(10) == 1) || (digitalRead(11) == 1)){`
- `Serial.println('!');`
- `//serial.println(analogRead(A0);`
- `}`

- else{
- // send the value of analog input 0:
- Serial.println(analogRead(A0));
- }
- //Wait for a bit to keep serial data from saturating
- delay(1);
- Serial.println(analogRead(A0));
- }

SOURCE OF PULSE-OXIMETER:

- `#include <Wire.h>`
- `#include "MAX30105.h"`
- `MAX30105 particleSensor;`
- `void setup() {`
- `Serial.begin(9600);`
- `// Initialize sensor`
- `if (particleSensor.begin() == false) {`
- `Serial.println("MAX30102 was not found. Please check`
`wiring/power.");`
- `while (1);`
- `}`
- `particleSensor.setup(); //Configure sensor. Use 6.4mA for LED drive`
- `}`

- void loop() {
- Serial.print(" R[");
- Serial.print(particleSensor.getRed());
- Serial.print("] IR[");
- Serial.print(particleSensor.getIR());
- Serial.println("]");
- }

APPLICATIONS AND ADVANTAGES OF ECG:

➤ **APPLICATIONS:**

- 1.fitness and activity heart rate monitors
- 2.portable ECG remote health monitors
- 3.gaming peripherals
- 4.biopotential signal acquisition

➤ **ADVANTAGES:**

- ECG helps to prevent heart attacks by analyzing heart parameters at the initial stage.
- ECG is used to detect the cardiac conditions of the patients after surgical or any other operation and after application of anesthesia.
- ECG test is quick, painless and safe.
- ECG test is cheap in cost.

APPLICATION AND ADVANTAGES OF PULSE-OXIMETER:

➤ **APPLICATION:**

- The pulse oximeter has already found a number of clinical applications outside of the operating room, such as monitoring during patient transport, respiratory monitoring during narcotic administration, and evaluation of home-oxygen therapy.

➤ **ADVANTAGES:**

- monitoring oxygen saturation over time.
- alerting to dangerously low oxygen levels, particularly in newborns.
- offering peace of mind to people with chronic respiratory or cardiovascular conditions.
- assessing the need for supplemental oxygen.

ADVANCEMENTS OF ECG AND PULSE-OXIMETER:

➤ **ECG:**


- Data can be uploaded or directly sent to the doctor by using IOT technology
- By using more electrodes measurements can be more accurate

➤ **PULSE OXIMETER:**

- Obtaining accurate pulse-oximetry readings in motion and low perfusion states, the development of a central sensor probe, miniaturization, and wireless technology are the future focus points for the ongoing development of pulse oximetry

EXPECTED RESULT


- **RESULT OF ECG:** If the test is normal, it should show that your heart is beating at an even rate of 60 to 100 beats per minute.
- Many different heart conditions can show up on an ECG, including a fast, slow, or abnormal heart rhythm, a heart defect, coronary artery disease, heart valve disease, or an enlarged heart.



EXPECTED RESULT

- **RESULT OF PULSE-OXIMETER:** A normal level of oxygen is usually 95% or higher. Some people with chronic lung disease or sleep apnea can have normal levels around 90%.
- The “SpO2” reading on a pulse oximeter shows the percentage of oxygen in someone's blood. If your home SpO2 reading is lower than 95%, call your health care provider.

Example of a pulse oximeter



CONCLUSION

- The 12-lead surface ECG can indicate pathological changes even before structural changes in the heart can be diagnosed by other methods. The recording of an ECG was of great value for several past generations of cardiologists and continues to provide vital information.
- pulse oximetry is a simple and non-invasive way to measure blood oxygen levels and heart rate. These measurements can be used to help monitor general health and quickly assess people with lung and heart disorders.

A central graphic featuring a teal watercolor splash or cloud-like shape on a light gray background. The splash has various shades of teal and blue, with a soft, painterly texture. The words "Thank You" are written in a black, cursive script font across the middle of this splash.

Thank You